

Medical Plaster Saw

S-750



Performance & Application

The Electric plaster saw S-750 is designed to remove all kinds of orthopedic plaster bandages etc. It is perfectly suited for the efficient removal of plaster and plastic materials. It boasts a robust construction, providing powerful and rapid cutting while ensuring skin safety. With its safe operation, low noise levels, and portability, this saw offers a superior experience. It effectively cuts all types of bandages, thanks to its strong motor and smooth, quiet motion. The saw includes automatic overload protection for added safety and durability during extended working hours. Its ergonomic design enables comfortable handling and easy maintenance.

- Plaster Saw
- Allen Key
- Round 65mm dia. Blade
- Round 42mm dia. Blade
- Universal Plaster Cutting Blade
- Storage/Carry Case



Packing



Aluminum case packing

Main Technical Data

Saw handpiece	Output voltage	80W	Non-loading speed	0-13000cpm
	Temperature rise	≤25° C	Non-loading noise	≤65dB (A)
	Output power	AC220V/50Hz or 110V/60Hz	Classification	Short-time duty