

Micro Medical Electric Spine Drill

S-720




Features

- Refined design
- Low noise, high speed
- Imported motor, autoclavable
- Long service life

Performance & Application


- **Straight handle:** different length, applied in neurosurgery, spine surgery, and microsurgery.
- **Angled handle:** different length, applied in neurosurgery, spine surgery, and microsurgery.



Packing & Sterilization Case



Aluminum case packing



Sterilization case
for power tool

Main Technical Data

Hand piece	output power	≥50W	Non-loading speed	0-40000 RPM
	Temperature rise	≤25° C	Non-loading noise	≤50dB (A)
Charger	Input Voltage	AC 220V/50Hz	Axial run out	≤0.03mm
			Maximum torque	0.5 N*m